


# ESTELITE OMEGA



# Concept : Simple & Excellent

- The best shade structure for layering technique
  - ✓ Optimization of value, chroma and transparency
  - ✓ High opalescence (Trans Shade)
- Simple system taking advantage of blending effect
- Excellent aesthetics


## Collaborated with:

Dr. Newton Fahl, Jr.



# Shade Structure

|                                | Shade Name  | Abb.  |
|--------------------------------|-------------|-------|
| Body Enamel                    | Enamel A1   | E A1  |
|                                | Enamel A2   | E A2  |
|                                | Enamel A3   | E A3  |
|                                | Enamel B1   | E B1  |
| Effect Translucent Enamel      | Trans       | Trans |
| Effect Semi-Translucent Enamel | Milky White | MW    |
| Dentin                         | Dentin A1   | D A1  |
|                                | Dentin A2   | D A2  |
|                                | Dentin A3   | D A3  |
| Bleach                         | Bleach 1    | BL 1  |
|                                | Bleach 2    | BL 2  |



# Delivery Form

## PLTs

- PLT DELUXE KIT
- PLT BASIC KIT
- Refill



## Syringes

- SYRINGE DELUXE KIT
- SYRINGE BASIC KIT
- Refill



# SYRINGE DELUXE KIT



Brush

Clinical Cases (DVD)  
Clinical Guide (CD)  
Restoration Techniques (Card)

Custom Shade Guide

- 11 syringes  
EA1, EA2, EA3, E B1  
DA1, DA2, DA3  
MW, Trans, BL 1, BL 2
- Custom Shade Guide x1
- Brush x1
- Clinical Cases (DVD) x1
- Clinical Guide (CD) x1
- Restorative Techniques (Card) x1
- Instruction x1

# SYRINGE BASIC KIT




Brush

Custom Shade Guide

Clinical Cases (DVD)  
Clinical Guide (CD)  
Restoration Techniques (Card)

- 6 syringes  
EA1, EA2, EB1  
DA2  
MW, Trans
- Custom Shade Guide x1
- Brush x1
- Clinical Cases (DVD) x1
- Clinical Guide (CD) x1
- Restorative Techniques (Card) x1
- Instruction x1

# PLT DELUXE KIT



Custom Shade Guide

Clinical Cases (DVD)  
Clinical Guide (CD)  
Restoration Techniques (Card)

- 11shade, 100 PLT  
Enamel A1,A2,A3,B1 (each 10PLT)  
Dentin A1,A2,A3 (each 10PLT)  
Milky-White (5PLT)  
Trans (5PLT)  
Bleach1,2 (each 10PLT)
- Custom Shade Guide x1
- Brush x1
- Clinical Cases (DVD) x1
- Clinical Guide (CD) x1
- Restorative Techniques (Card) x 1
- Instruction x1

# PLT BASIC KIT



Brush

Custom Shade Guide

Clinical Cases (DVD)  
Clinical Guide (CD)  
Restoration Techniques (Card)

- 6shade, 50 PLT  
Enamel A1,A2,B1 (each 10PLT)  
Dentin A2 (10PLT)  
Milky-White (5PLT)  
Trans (5PLT)
- Custom Shade Guide x1
- Brush x1
- Clinical Cases (DVD) x1
- Clinical Guide (CD) x1
- Restorative Techniques (Card) x 1
- Instruction x1

# Syringe Refill



- 4g syringe x 1
- Instruction x 1


# PLT Refill



- 0.25g PLT x 20
- Instruction x 1

# Clinical Guide (CD)

Step-by-step description by Dr. Newton Fahl, Jr.



**CLASS IV RESTORATIONS**

Class IV restorations can be performed with ESTELITE OMEGA according to three layering approaches depending on the type of translucency pattern to be emulated. Three shades are always present in all three: a thin layer of Milky-White is used to establish a lingual frame, Dentin is used for opacity, and Body Enamel generates outer color and a certain degree of translucency. When used internally at the incisal one third, Trans Enamel promotes deeper translucency and opalescence. For further translucent effects, Milky-White can be applied incisally over the Trans Enamel.

**4-layer approach – step-by-step**

Fig. 1 – A midline Diastema of moderate size (1.8 mm) – labial view.  
Fig. 2 – A midline Diastema of moderate size (1.8 mm) – facial view.  
Fig. 3 – Technical approach is limited to removal of the apical enamel layer through prophylaxis or sandblasting.

**3-layer approach** (minimal translucency and dentin mamelon effects at incisal third)

Fig. 4 – An Estelite shade is selected to match the class of the cervical and middle thirds.  
Fig. 5 – Milky-white is selected to be used both as a lingual shell and a final varnish over the incisal one third.  
Fig. 6 – Enamel stain is selected to impart greater depth and opalescence where needed incrementally.

**4-layer approach** (moderate translucency and dentin mamelon effects at incisal third)

Fig. 7 – A modified rubber-dam technique and retraction cord are used for field control.  
Fig. 8 – The enamel surface is etched for 15-30 seconds with phosphoric acid.  
Fig. 9 – Bond force is applied and gently rubbed against the enamel surface for 20 seconds.

**5-layer approach** (maximum translucency and dentin mamelon effects at incisal third)

Contents

Introduction  
Class IV  
Direct Veneer  
Diastema Closure  
Class V  
Class I & II

# Restorative Techniques (Card)

## Extract from Clinical Guide

**ESTELITE OMEGA Restorative Techniques**  
by Dr. Newton Fahl Jr.

**CLASS IV RESTORATIONS**  
Class IV restorations can be performed with ESTELITE OMEGA according to three layering approaches depending on the type of translucency pattern to be simulated. Three shades are always used: a thin layer of Milky-White is used to establish a gingival trans. Dentin is used for opacity, and Body Enamel generates color and a certain degree of translucency. When used internally at the incisal one third, Trans Enamel provides deeper translucency and opalescence. For further translucent effects, Milky-White can be applied incisally over the TransEnamel.

**DAIStEMA CLOSURE**  
Diastema closures can be performed according to two restorative approaches. The 2-layer technique is indicated for cases where there is little or no translucent effect at the proximal aspect of the incisal one third and requires a Milky-White Enamel shell and a Body Enamel as a fill layer. The 4-layer approach is indicated when opalescence associated with a thin translucency is desired at the incisal one third such effect, a Trans Enamel is used and is covered by a Milky-White Enamel further promote light transmission and

**CLASS V RESTORATIONS**  
Class V defects can be restored via a monochromatic or bichromatic approach. Most Class V restorations require the use of higher opacity/higher chroma shades because of the opacity/color characteristics.


**ESTELITE OMEGA Restorative Techniques**  
by Dr. Newton Fahl Jr.

**DIRECT VENEER RESTORATIONS**  
Direct veneers can be classified in three types according to the need for value change and for incisal edge modification. Type 1 is a basic approach for cases requiring no incisal edge modification and when the dentition presents with a natural basic color but the form needs improvement. Type 1 involves one layer of Body Enamel and one layer of Milky-White Enamel. If a higher chroma is desired in specific areas, such as the cervical one third, ESTELITE COLOR (Ochre, Yellow, Red) can be used. Types 2 and 3 entail the need for incisal edge length modification. Similarly to Type 1, Type 2 is meant for changing mostly the form of a tooth that is not discolored, only requiring the need of an ESTELITE OMEGA Dentin for added opacity and Body Enamel for the basic hue, and Trans, and Milky-White Enamel for determining varying degrees of translucency effects. Type 3 also requires modifying the incisal edge and is for more complex cases needing masking of tooth discoloration. ESTELITE COLOR (Low Chroma Opaque, Medium Chroma Opaque, High Chroma Opaque) is used to block out the dark tooth substrate and to elevate the value to a more luminous level. The remaining layers are much like those of Type 2 and are used to provide beautiful color nuances with dentin and enamel shades, including a bleach shade to create white characterizations. A variation of Type 3 is also possible, whereby no incisal edge modification is realized. In the diagrams below, the colored bullets depict the shades used and the sequence according to which each layer is applied.

**CLASS I and II RESTORATIONS**  
Class I and II restorations can be performed with ESTELITE OMEGA according to two basic layering approaches: Monochromatic and Bichromatic. Monochromatic restorations require either an A1 or A2 shade of Body Enamel for a beautiful color blend in. Bichromatic restorations require a higher chroma Dentin (A3) and a Milky-White Enamel for a restoration of superior aesthetics. ESTELITE COLOR tints can optionally be used to further import polychromaticity.


**Monochromatic Approach**  
(1 Dentin + Milky-White Enamel, 1 Dentin)

**Bichromatic Approach**  
(1 Dentin + Milky-White Enamel, 1 Dentin)



Contents  
Introduction  
Class IV  
Direct Veneer  
Diastema Closure  
Class V  
Class I & II


# Custom Shade Guide




Mold x 1




Tips x 15




Holder x 1




Set a tip to mold



Fill resin paste to mold



Light cure



# Brush



Flat Brush

Thin Brush

# ESTELITE COLOR

13 Shades

Clear  
White  
Blue  
Yellow  
Gray  
Ochre  
Dark Brown  
Red  
Lavender

Pink Opaque  
High Chroma Opaque  
Medium Chroma Opaque  
Low Chroma Opaque

| shade  | 0.1mm |            |
|--------|-------|------------|
| Clear  |       | Ochre      |
| White  |       | Dark Brown |
| Blue   |       | Red        |
| Yellow |       | Lavender   |
| Gray   |       |            |

# ESTELITE COLOR

## KIT

0.9g syringe x 13shades  
Syringe Tip & Cap x 15  
Disposable Dish x 5  
Light Shield Cover x 1  
Instruction

## Refill


0.9g syringe x 1  
Instruction

# Clinical Cases

Clinical Case

# Class IV Restoration

Dr. Newton Fahl Jr.



Clinical Case

# Direct Veneer


Dr. Newton Fahl Jr.



Clinical Case

# Diastema Closure

Dr. Newton Fahl Jr.



Clinical Case

# Class V Restoration

Dr. Newton Fahl Jr.

