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PF

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PHS DO BRASIL LTDA

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Versão: 09

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PF

INSTRUCTIONS FOR USE FOR PROFESSIONALS

PHS DO BRASIL LTDA

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: 8.11992.7 (X401M5619X7L)

Notification at ANVISA number: 8.11992.70006

EU

REP

CINTERQUAL – Commerce Solutions

International Lda; Avenida Defensores de Chaves, Nº 4

Idea Spaces Office – 1000-117

Lisboa – Portugal +351 967 417 654

Product Commercial Name: POTENZA BIANCO CARBAMIDE PEROXIDE PF 10%, 16%, 22%.

Alternative Commercial Name: POTENZA BIANCO PF DENTAL WHITENING FOR HOME USE.

Technical Product Name: DENTAL WHITENING AGENTS.

READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.

KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.

Technical Support

0800-000-3780

www.phsgroup.com.br

BASIC COMPOSITION:

Carbamide peroxide, carbomer, potassium nitrate, sodium fluoride, humectant (glycerol), neutralizing agents and ultrapure water.

PRODUCT DESCRIPTION:

POTENZA BIANCO PF is a dentist-supervised at-home dental whitening agent available in 10%, 16%, and 22% Carbamide Peroxide concentrations. The product is supplied as a colorless thixotropic gel packaged in ready-to-use syringes.

Its formulation contains Potassium Nitrate and Sodium Fluoride, agents widely used to help reduce tooth sensitivity associated with whitening procedures. In addition, it has a high water content and contains Glycerol, contributing to the maintenance of gel hydration throughout the treatment.

PRODUCT INDICATION:

POTENZA BIANCO PF is indicated for dentist-supervised at-home dental whitening. It may be used on vital and non-vital teeth (external application only) presenting discoloration. The product is also indicated as an adjunctive treatment for mild cases of fluorosis staining or staining caused by tetracycline and its derivatives.

PRECAUTIONS:

- If the patient reports pre-existing tooth sensitivity during daily activities (e.g., consumption of cold beverages, exposure to cold air, or pressure applied to the teeth) and no apparent clinical cause is identified, it is recommended to perform a desensitizing protocol with Potenza Esente 0.2% and Potenza Esente GLU HEMA before whitening treatment;
- In cases of significant gingival recession with root exposure, it is recommended to restore or protect the exposed area with an appropriate material, according to the dentist's clinical judgment;
- Dental surfaces that will receive the whitening gel must be clean and free from dental plaque and calculus. When necessary, prophylaxis and/or scaling procedures should preferably be performed at least 5 days before starting the whitening treatment;
- It is essential that the whitening gel has minimal contact with soft tissues, as it may cause local irritation of varying intensity depending on exposure time, product concentration, and individual tissue characteristics. Instruct the patient to apply a small amount of whitening gel (approximately half a grain of rice) into the corresponding space for each tooth in the tray, avoiding excessive overflow beyond the tray margins;
- If the patient experiences a high degree of tooth sensitivity during treatment, the clinical protocol should be reassessed. When applicable, reduce the whitening gel concentration to 10% and instruct the patient to use Potenza Esente 0.2% before and after each whitening application for 10 minutes.

IMPORTANT:

- Whitening treatment must be carried out exclusively under the supervision of a dentist;
- The safety and effectiveness of the treatment depend on proper indication and use of the product according to professional guidance and these Instructions for Use;
- Whitening effectiveness may vary among patients due to individual factors such as the type and intensity of dental stains, dental tissue permeability, dentin characteristics, and other individual clinical conditions;
- In certain cases, the dentist may combine the treatment with in-office whitening procedures according to clinical evaluation.

CONTRAINDICATIONS:

Before starting whitening treatment, perform a thorough patient assessment, including medical and dental history, habits, systemic conditions, and possible risk factors. During the clinical examination, identify any conditions that require prior treatment or contraindicate whitening procedures.

A) TEMPORARY CONTRAINDICATIONS – Must be resolved before treatment:

- Presence of dental caries;
- Defective or leaking restorations;
- Dental cracks or fractures;
- Periodontal diseases, including gingivitis and periodontitis;
- Presence of dental plaque and/or calculus;
- Teeth with questionable pulpal vitality;
- Untreated endodontic conditions;
- Extensive gingival recession with significant root exposure;
- Severe tooth sensitivity without prior clinical investigation.

B) ABSOLUTE CONTRAINDICATIONS – Do not perform treatment:

- Patients with a history of hypersensitivity or allergy to any component of the formulation;
- Teeth presenting amelogenesis imperfecta or dentinogenesis imperfecta;
- Teeth with incomplete root development;
- Teeth with severe mobility or severe periodontal compromise;
- Pregnant or breastfeeding women;
- Individuals under 15 years of age;
- Patients with severe uncontrolled systemic diseases or those for whom medical or dental evaluation contraindicates treatment.

POSSIBLE SIDE EFFECTS:

- It is common for the patient to present a certain degree of hypersensitivity during or after the application of the whitening gel. This is due to the product's own mechanism of action. The level of

discomfort varies greatly from patient to patient. The desensitizing agent Potenza Esente 0.2% should be associated with the treatment for greater patient comfort. In extreme cases, it is suggested to stop the treatment;

- The bleaching can highlight pre-existing areas of hypocalcification (hypoplasia) due to dehydration or even more intense bleaching of these areas. In most cases, these spots disappear after rehydration or remineralization of the region, which occurs within a few hours or days after using the bleaching agent. Whenever you identify the stains, suspend the treatment until they have disappeared. If the stains are persistent, discuss the continuation or not of the treatment with the patient, informing him of any need for restorative intervention;
- Some patients may have non-uniform whitening in different thirds of the teeth. This is due to a physiological issue and will hardly be changed, because in the cervical region of the teeth there is a greater thickness of dentin, which is naturally more stained than the predominant enamel in the other thirds.

WARNING:

- Observe the expiration date indicated on the packaging. Do not use the product after expiration;
- Keep out of reach of children;
- Product intended exclusively for dental use;
- Before use, inspect the packaging integrity and product appearance. If any color change, consistency alteration, leakage, or other abnormality is observed, discontinue use and contact the manufacturer;
- The product must be used exclusively according to the dentist's instructions and these Instructions for Use;
- Avoid prolonged contact of the gel with gingiva, oral mucosa, lips, or skin, as temporary local irritation may occur;
- In case of accidental eye contact, rinse immediately with plenty of water and seek medical attention;
- Do not ingest the product. In case of accidental ingestion, seek medical attention;
- If persistent irritation, excessive sensitivity, or any unexpected adverse reaction occurs, discontinue use and consult a dentist;
- The manufacturer is not responsible for damages resulting from improper use, use contrary to these Instructions for Use, or incorrect handling of the product;
- It is the responsibility of the professional to evaluate the clinical indication and suitability of the product for each patient before use;
- Do not reuse the empty packaging for any purpose;
- Dispose of the product and its packaging according to local regulations.

INSTRUCTIONS FOR USE:

Before using the product, carefully read all information contained in this document. Perform a complete medical and dental assessment and resolve any conditions that may contraindicate or interfere with whitening treatment (see Contraindications and Precautions sections).

SUGGESTION OF USAGE TIME ACCORDING TO GEL CONCENTRATION
CARBAMIDE PEROXIDE / DAILY USE TIME

10%	3 TO 4 HOURS
16%	2 TO 3 HOURS
22%	1 HOUR

Attention: The dentist must monitor the patient throughout the entire at-home treatment. Periodic follow-up appointments are recommended to evaluate whitening progress and possible adverse effects.

1st Step: Select and record the patient's initial tooth shade using an appropriate shade guide. Photographic documentation is recommended for monitoring clinical progress.

2nd Step: Obtain the patient's dental models through an appropriate impression-taking procedure. Trim the models into a horseshoe shape, removing the palate reproduction from the upper model and the tongue reproduction from the lower model.

3rd Step: Position the model in the thermoforming machine. Install a 1.25 mm EVA sheet (Potenza Stampo) and start the heating cycle. When sufficient softening of the sheet occurs (approximately 2 cm of sagging), activate the vacuum system to form the tray over the model.

4th Step: Allow the sheet to cool until it reaches room temperature. Remove the tray from the model and trim it as necessary. It is recommended that the tray borders remain approximately 1 mm short of the gingival margin. Remove any edges or sharp points that may cause discomfort to the soft tissues.

5th Step: Instruct the patient to apply the whitening gel directly from the syringe using the application tip. The gel should be placed in the area corresponding to each tooth to be whitened, typically from second premolar to second premolar, using an amount approximately equivalent to half a grain of rice per tooth, positioned in the center of the facial surface area of the tray.

6th Step: Instruct the patient to properly seat the tray onto the dental arch. Any excess gel that extrudes beyond the tray margins should be removed immediately using clean fingers or a moistened cotton swab.

7th Step: During the tray wearing period, the patient must not eat, drink, or smoke.

8th Step: After the recommended wearing time, remove the tray and thoroughly rinse the mouth with water. Perform normal oral hygiene procedures.

9th Step: Wash the tray with running water at room temperature. After removing all gel and saliva residues, dry it with a paper towel and store it in an appropriate case, protected from moisture, excessive heat, and direct sunlight.

NOTE:

If the patient reports frequent swallowing of whitening gel while wearing the tray or presents signs of gingival irritation, consider reducing the tray border extension and instruct the patient to apply a smaller amount of gel to each tooth.

GENERAL SAFETY INSTRUCTIONS:

- The carbamide peroxide contained in the syringe may be under pressure and/or may have leaked due to adverse storage conditions. Before opening the package, protect your hands and eyes and inspect the package for any signs of leakage;
- Avoid contact of the product with skin, mucous membranes, or clothing;
- Instruct the patient not to swallow the product. While wearing the tray loaded with whitening gel, the patient should not eat food or drink water;
- In case of accidental contact of the gel with the skin or oral mucosa, rinse the affected area thoroughly with water;
- In case of accidental contact of the whitening gel with the eyes, rinse immediately with plenty of water and seek medical attention, taking these Instructions for Use with you;
- In the event of allergic reactions, treatment must be discontinued immediately. Refer the patient for medical evaluation and provide these Instructions for Use;
- Leftover or residual whitening gel must not be reused. Dispose of the product according to local regulations.

CONSERVATION AND STORAGE:

- Keep the product in its original packaging;
- Protect from direct sunlight;
- Store the product at temperatures from 0 °C to 30 °C (32 °F to 86 °F).

TRANSPORT CONDITIONS:

- The product can be transported at temperatures between -5 °C and 50 °C (23 °F to 122 °F);
- Protect from direct sunlight.

MANUFACTURING DATE: SEE PACKAGING.

EXPIRATION DATE: 2 YEARS AFTER THE MANUFACTURE DATE.

SALES UNDER PRESCRIPTION (ACCORDING TO RDC N° 6 OF 2015).

MAKE SURE THE INSTRUCTION FOR USE MATCHES THE PRODUCT TO BE USED.

TO ACCESS A PRINTED VERSION OF THE INSTRUCTION FOR USE, CONTACT THE COMPANY'S TECHNICAL SUPPORT.

KEEP THE PACKAGE CLOSED WHEN NOT IN USE.

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