# Potenza BLOXY

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## **PHS DO BRASIL LTDA**

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# Potenza BLOXY

INSTRUCTIONS FOR USE FOR PROFESSIONALS

## **PHS DO BRASIL LTDA**

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: X401M5619X7L (8.11992.7)

Notification at ANVISA number: 8.11992.70030

**Product Commercial Name:** POTENZA BLOXY.

**Alternative Commercial Name:** OXYGEN BLOCKING GEL.

**Technical Product Name:** DENTAL AUXILIARY MATERIALS.

**READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.**

**KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.**

**BASIC COMPOSITION:**

Glycerin, thickener, preservative, neutralizer, chelating agent and aroma.

**PRODUCT DESCRIPTION:**

The operating principle involves creating a physical, impermeable barrier between dental materials, such as composite resins and resin cements, and oxygen from the surrounding environment. This is essential to avoid inhibiting the polymerization of these materials by oxygen, which could compromise the quality and effectiveness of dental procedures. Glycerin, as a key component of the gel, plays an important role in this oxygen-locking process. Glycerin is a viscous liquid that has lubricating properties and is capable of forming a protective layer on the surface of dental materials. The application of the glycerin-based gel forms a physical barrier between the material and oxygen, preventing oxygen from reaching the material and interfering with the polymerization reaction. Precisely due to the gel's characteristics, it is also capable of lubricating dental materials such as gingival retraction threads, in a safe and gentle way.

**PRODUCT INDICATION:**

POTENZA BLOXY has lubricating properties and is resistant to oxygen penetration, being formulated with a glycerin base for the purpose of establishing a protective barrier. This barrier prevents contact between composite resins and resin cements with oxygen, ensuring greater polymerization effectiveness. Furthermore, this gel plays a role in lubricating the gingival retraction threads, providing a delicate and efficient application.

**PRECAUTIONS:**

- Sensitivity Test: Before applying the gel to a patient, perform a sensitivity test to check for allergic reactions or irritations. Apply a small amount of the gel to a discreet area and observe possible reactions for a short period;
- Avoid direct contact of the gel with the eyes and skin. Use adequate protection, such as gloves and protective glasses, during application;
- Apply the appropriate amount of gel to cover the desired areas. Avoid excesses that could make the procedure difficult or interfere with polymerization;
- Remove the oxygen blocking gel completely and efficiently after the polymerization step. Ensure that no residue remains in the treated area;
- Observe the patient during and after the procedure to detect any adverse reactions, such as irritation, allergies or discomfort;
- Check the compatibility of the gel with other dental materials, such as resins and cements, to avoid interference with adhesion or physical properties;
- Keep the oxygen blocking gel in a cool, dry place, away from direct heat sources and direct sunlight.

**CONTRAINDICATIONS:**

- Oxygen blocking gel may not be suitable for young children, especially those who have difficulty following instructions or who may accidentally swallow the product;
- If the patient has had previous adverse reactions to similar dental products, it is prudent to avoid the use of oxygen blocking gel to minimize the risk of recurrence;
- Do not use the material after the expiration date indicated on the packaging.

**POSSIBLE SIDE EFFECTS:**

- Some patients may experience irritation or sensitivity in the area where the gel was applied. This may manifest as redness, swelling, burning or discomfort;
- Although rare, some people may have allergic reactions to ingredients in oxygen blocking gel. This can lead to itching, rashes, swelling, or other allergic symptoms;
- If the gel is not applied correctly or removed properly after photoactivation, interference with the polymerization of dental materials may occur, affecting the quality of the treatment.

**WARNING:**

- Dental Oxygen Blocking Gel is intended solely for use by trained dental health professionals. It should not be used by patients without adequate guidance or supervision;
- Do not ingest the product;
- Avoid direct contact with eyes. In case of accidental contact, wash eyes with plenty of water and seek medical assistance if irritation persists. Also, avoid prolonged contact with the skin and wash your hands after application;
- Check the compatibility of the gel with other dental materials, such as resins, cements and devices. Some interactions may affect adhesion or mechanical properties;
- Dispose of any excess oxygen blocking gel in accordance with local regulations. Avoid environmental contamination and exposure to third parties;
- Respect the expiration date. Never use the product if it has expired;
- Keep out of the reach of children;
- This product was manufactured for dental use only and must be handled in accordance with the instructions for use;
- It is the professional's responsibility to check, before using the product, whether this product is compatible and suitable for the desired use;
- Do not reuse empty packaging for other purposes;
- PRODUCT FOR EXCLUSIVE USE IN DENTAL OFFICES;
- ONLY FOR USE BY DENTAL PROFESSIONALS;
- KEEP THE PACKAGE CLOSED WHEN NOT IN USE;
- IMPORTANT: Suspend use of the product in case of allergic reactions.

## **INSTRUCTIONS FOR USE:**

### **Cementation of indirect parts:**

**1st STEP:** With the resin cement already applied inside the piece, take it to the tooth and remove excess cement that leaks;

**2nd STEP:** Apply a layer of Potenza Bloxy to the entire accessible margin of the piece and light-polymerize the entire region for the time recommended by the cementing agent manufacturer. If you used a self-polymerizing cement, keep Potenza Bloxy applied throughout the product's chemical curing time;

**3rd STEP:** Remove Potenza Bloxy with air/water spray for 30 seconds or until completely removed. Then proceed to finish and polish the margins as necessary.

### **Direct Restores:**

**1st STEP:** After completing the sculpting of the last increment of the restoration and before light-curing it, apply a thin layer of Potenza Bloxy gently over the entire accessible surface of the restoration and its visible margins. Be careful not to change the shape of the restoration during application;

**2nd STEP:** With the surface entirely covered, light cure for the time recommended by the composite manufacturer;

**3rd STEP:** Remove Potenza Bloxy with air/water spray for 30 seconds or until completely removed. Then proceed to finish and polish the restoration as necessary.

## **CONSERVATION AND STORAGE:**

- Keep the product in the original packaging;
- Protect from direct sunlight;
- Store the product at temperatures from 0 °C to 30 °C.

## **TRANSPORT CONDITIONS:**

- The product can be transported at temperatures between -5 °C and 50 °C.

**MANUFACTURING DATE: SEE PACKAGING.**

**SHELF LIFE: 2 YEARS AFTER THE MANUFACTURING DATE.**

**PRODUCT OF EXCLUSIVE USE IN DENTAL OFFICE. ONLY FOR USE OF DENTAL PROFESSIONALS.**

**MAKE SURE THAT THE INSTRUCTION FOR USE IS IN ACCORDANCE WITH THE PRODUCT TO BE USED.**

**TO ACCESS THE PRINTED VERSION OF THE INSTRUCTION FOR USE, CONTACT THE COMPANY'S TECHNICAL SUPPORT.**

**KEEP THE PACKAGING CLOSED WHEN NOT IN USE.**

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