TOKUYAMA BRUSH Nº24

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PHS DO BRASIL LTDA

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Versão: 05

TOKUYAMA BRUSH Nº24

INSTRUCTIONS FOR USE FOR PROFESSIONALS

PHS DO BRASIL LTDA

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: X401M5619X7L (8.11992.7)

Notification at ANVISA number: 8.01175.80435

Product Commercial Name: BRUSH TOKUYAMA Nº 24.

Technical Product Name: BRUSH FOR MODELING LIGHT-CURED COMPOSITE RESIN.

READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.

KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.

Technical Support

0800-000-3780

www.phsgroup.com.br

PRODUCT'S DESCRIPTION:

Nylon bristles, metal ferrule/ferrule and painted/varnished wooden handle with print.

DESCRIPTION AND INDICATIONS:

Tokuyama Brush n° 24 is a brush with two distinct active ends, both synthetic (nylon), one with a flat morphology and the other tapered. The brush is recommended for dental use during the modeling and sculpting stage of composite resins in direct restorations. The flat end is recommended for blending and smoothing layers, while the tapered end is for texturing and detailing.

PRECAUTIONS AND CONTRAINDICATIONS:

- Keep away from heat and/or sunlight;
- Keep out of reach of children;
- Do not use on patients allergic to any of the product's components;
- In cases of allergy, refer the patient to medical care;
- Non-injectable product;
- Do not ingest the product;
- Do not use for purposes other than those described in this manual.

POSSIBLE SIDE EFFECTS:

- No known side effects if used according to instructions for use.

WARNING:

- For dental use only;
- This product must be handled in accordance with these instructions for use;
- It is the professional's responsibility to check before using the product whether it is compatible and suitable for the intended use;
- Before using, observe the appearance of the product; if you notice any anomaly, contact the manufacturer immediately;
- The manufacturer is not responsible for damages caused by inappropriate use or incorrect handling;
- To dispose of the product, follow local legislation.

INSTRUCTIONS FOR USE:

Before and after each use of the product, follow the steps as recommended in "Conservation and Storage" so that the brush is clean and decontaminated.

1st STEP: Using spatulas, apply an increase in the composite resin to the dental substrate, modeling its basic shape;

2nd STEP: Using the flat end of the brush, brush the resin increment against the tooth surface until the layer blends with the tooth or with the pre-existing layers of resin. This makes the restoration surface smoother and more homogeneous, facilitating the finishing and polishing step later;

3rd STEP: With the tapered end of the brush, it is possible to sculpt anatomical details of the restoration (grooves, angles, edges) as needed. This end also allows you to apply dyes or the like for characterization.

CONSERVATION AND STORAGE:

To clean and disinfect the instrument, follow these steps:

- Gently rub a gauze pad soaked in alcohol 70 in gentle movements in the direction of the bristles. Repeat until resin/adhesive residue is removed;
- To remove alcohol residue, repeat the maneuver but with gauze soaked in water;
- For disinfection, submerge the active tip of the brush in a 2% chlorhexidine solution for 5 minutes. Then, use deionized water to remove residues from the solution;
- To remove excess water from the bristles, brush gently on sterile gauze. Adapt the protector (plastic cylinder) to the active tip and store it in a dry environment free from direct sunlight. Alternatively, the bristles can receive glycerin gel before storing the instrument, remembering to remove it before the next use.

STORAGE CONDITIONS:

- Store the brush in a dry place, protected from direct sunlight, at room temperature.

TRANSPORT CONDITIONS:

- The product can be transported at temperatures between -5 °C and 50 °C;
- Protect from direct sunlight.

MANUFACTURING DATE: SEE PACKAGING.

EXPIRATION DATE: UNDETERMINED.

MAKE SURE THE INSTRUCTION FOR USE MATCHES THE PRODUCT TO BE USED.

TO ACCESS A PRINTED VERSION OF THE INSTRUCTION FOR USE, CONTACT THE COMPANY'S TECHNICAL SUPPORT.

KEEP THE PACKAGE CLOSED WHEN NOT IN USE.

PRODUCT OF EXCLUSIVE USE IN DENTAL OFFICE. ONLY FOR USE OF DENTAL PROFESSIONALS.

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