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PHS DO BRASIL LTDA

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Código: IFU-039

Versão: 00

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INSTRUCTIONS FOR PROFESSIONAL USE

PHS DO BRASIL LTDA

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: 8.11992.7 (X401M5619X7L)

Notification at ANVISA number: 81199270035

Product Commercial Name: Potenza FineX Cut

Alternative Commercial Name: Potenza FineX Cut Interproximal Saw

Technical Product Name: Interproximal Saw

READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.

KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.

Technical Support

0800-000-3780

www.phsgroup.com.br

BASIC COMPOSITION:

- 54 mm interproximal saws made of AISI 304 stainless steel, provided with autoclavable grip buttons at both ends.
- 54 mm interproximal saws made of AISI 304 stainless steel with a titanium-based surface coating, provided with autoclavable grip buttons at both ends.
- 150 mm interproximal saws made of AISI 304 stainless steel.
- 150 mm interproximal saws made of AISI 304 stainless steel with a titanium-based surface coating.

PRODUCT DESCRIPTION:

The Potenza FineX Cut Interproximal Saw is a professional-use dental medical device consisting of a serrated blade manufactured from AISI 304 stainless steel. Only the 54 mm versions are provided with grip buttons at both ends to improve stability, comfort, and control during handling.

The product is marketed as a family of devices sharing the same design characteristics, materials, principle of operation, and risk profile, differing only in blade dimensions and surface finish. It is available in versions with or without a titanium coating. Titanium-coated models provide increased wear resistance and improved cutting performance.

The product family consists of 54 mm and 150 mm versions. Both are supplied non-sterile and must be cleaned and sterilized before use according to the instructions provided in this IFU.

PRODUCT' INDICATION:

The Potenza FineX Cut Interproximal Saws are intended for the manual removal of excess restorative materials and for interproximal finishing of dental restorative materials such as composite resins, cements, adhesives, and amalgam, providing refinement of interproximal surfaces without compromising the contact point between adjacent teeth.

The device is intended for professional use by dentists during restorative, prosthetic, and aesthetic dental procedures.

PRECAUTIONS:

- Visually inspect the interproximal saw before each use. Do not use the product if deformation, corrosion, fractures, excessive wear of the cutting teeth, or any other condition that may compromise its performance is observed.
- After each use, perform the cleaning, drying, and sterilization procedures as described in this Instructions for Use.
- During cleaning, do not use abrasive products, metal brushes, or incompatible chemical agents, as they may damage the instrument's surface.

- Before first use, verify that the package is intact and shows no signs of damage or contamination.
- Replace the product if significant loss of cutting efficiency, deformation, corrosion, fractures, or any condition affecting its physical integrity or performance is observed.
- After cleaning and sterilization, store the saw in a clean, dry location protected from moisture, corrosive agents, and direct sunlight, keeping it in an appropriate sterilization pouch until its next use.

CONTRAINDICATIONS:

The Potenza FineX Cut Interproximal Saw is contraindicated in the following situations:

- Patients with known hypersensitivity or allergic reactions to any of the product materials.
- Use for purposes other than those described in these Instructions for Use.
- Use in clinical situations where, in the dentist's judgment, the use of the saw may pose an additional risk to the patient or compromise the outcome of the procedure.
- Use of a product showing deformation, corrosion, fractures, or any condition that may compromise its physical integrity or performance.

The decision to use the product should always be based on the individual clinical assessment performed by the responsible dentist.

POSSIBLE SIDE EFFECTS:

During use of the Potenza FineX Cut Interproximal Saw, occasional effects related to the clinical procedure or mechanical contact between the instrument and oral tissues may occur, including:

- Patient discomfort or transient sensitivity during the procedure;
- Irritation or minor injuries to the gingiva or soft tissues caused by accidental contact with the saw;
- Minor gingival bleeding, especially in areas with pre-existing inflammation or during restorative procedures close to the gingival sulcus;
- Transient tooth sensitivity following finishing or adjustment of the restoration;
- Hypersensitivity reactions in patients allergic to the materials composing the product.

The occurrence of these effects may be minimized through appropriate clinical technique, correct use of the instrument, and compliance with the instructions contained in this IFU.

WARNING:

- For professional use by dentists or other legally qualified dental professionals only.

- Before use, verify during patient assessment that the patient has no known history of hypersensitivity to any of the product materials.
- Use the saw exclusively for the intended purpose described in these Instructions for Use.
- Before first use and between subsequent uses, the saw must be cleaned and sterilized as described in this IFU. Failure to properly clean and steam sterilize the product may result in cross-contamination.
- The use of a rubber dam or another appropriate isolation technique is recommended to reduce the risk of accidental ingestion or aspiration of the instrument or restorative material fragments.
- Wear appropriate personal protective equipment (PPE), including protective eyewear for both the practitioner and the patient, as well as a protective face mask during the procedure.
- Before starting the procedure, ensure that the operative field is adequately illuminated and that visibility is sufficient for safe execution of the technique.
- Maintain a firm grip on the instrument throughout the procedure and apply only light to moderate pressure. Avoid excessive force, twisting, or abrupt movements that may cause deformation, fracture of the saw, or injury to the patient.
- Avoid accidental contact between the active portion of the saw and the gingiva, oral mucosa, or other soft tissues, as this may result in local injury.
- Do not use the product if it shows deformation, corrosion, fractures, excessive wear of the cutting teeth, or any other condition that may compromise its physical integrity or performance.
- Immediately discontinue use if any instrument malfunction is observed during the procedure.

INSTRUCTIONS FOR USE:

1. Before first use and between subsequent uses, clean and sterilize the saw according to the instructions provided in this IFU.
2. Hold the interproximal saw firmly with both hands. When applicable (54 mm models), use the grip buttons to provide greater stability and control during the procedure.
3. Visually identify the cutting direction of the saw. The cutting teeth must be oriented away from the operator.
4. Identify the interproximal area to be treated and carefully position the saw in place. The saw should remain perpendicular to the long axis of the teeth.
5. Apply light to moderate pressure in the gingival direction while performing short, repetitive pulling strokes away from the operator until the excess restorative material (composite resin, cement, adhesive, amalgam, or other indicated materials) has been removed.
6. During the procedure, allow the holes in the body of the saw to collect cutting debris, thereby reducing accumulation in the active area.
7. After completing the cut, carefully remove the saw from the interproximal space.
8. Evaluate the quality of the interproximal finishing and the integrity of the contact point. If necessary, perform additional adjustments using appropriate instruments.

9. After completing the procedure, clean, dry, sterilize, and store the saw according to the instructions provided in this IFU.

CLEANING, DISINFECTION AND STERILIZATION

The Potenza FineX Cut Interproximal Saw is supplied non-sterile and must be cleaned and sterilized before first use and between subsequent uses.

Cleaning

- Immediately after use, remove all visible debris under running water.
- Wash the instrument using an enzymatic detergent according to the detergent manufacturer's instructions. If necessary, use a soft-bristled brush to remove adhered residues.
- Rinse thoroughly under running water until all detergent residues have been completely removed.
- Dry the instrument completely using a clean, lint-free cloth or oil-free, moisture-free compressed air.

Sterilization

- After drying, place the saw in a suitable steam sterilization pouch and seal it properly.
- Sterilize in a steam autoclave at 121 °C for 20 to 30 minutes or, alternatively, at 134 °C for 5 to 10 minutes, according to the validated sterilization protocol adopted by the healthcare facility.
- At the end of the sterilization cycle, allow the package to dry completely before removing it from the autoclave.
- Store the sterilized product in a clean, dry place, protected from moisture and direct sunlight, preserving the integrity of the package until use.

CONSERVATION AND STORAGE:

- Store the Potenza FineX Cut in a clean, dry, well-ventilated place, protected from moisture, dust, and corrosive chemical agents.
- Store the product at a temperature between 32°F and 104°F (0°C and 40°C).
- After cleaning, drying, and sterilization, keep the Potenza FineX Cut packaged in suitable steam sterilization packaging, preserving the integrity of the sterilization packaging and the sterility of the product until the time of use.

TRANSPORT CONDITIONS:

- The product must be transported in its original packaging in order to preserve its physical integrity and avoid deformation, mechanical impacts, perforations, contamination, and exposure to corrosive agents.

- During transport, the Potenza FineX Cut Interproximal Saw must be kept at a temperature between 0 °C and 60 °C, protected from excessive moisture and direct exposure to weather conditions.
- Stacking, when applicable, must respect the packaging resistance, avoiding excessive compression that may damage the product or its packaging.

DATE OF MANUFACTURE: SEE PACKAGING.

SHELF LIFE: SEE PACKAGING.

PRODUCT FOR EXCLUSIVE USE IN DENTAL OFFICES. ONLY FOR USE BY DENTAL PROFESSIONALS.


MAKE SURE THAT THE INSTRUCTIONS FOR USE CORRESPOND TO THE PRODUCT TO BE USED.

TO ACCESS THE PRINTED VERSION OF THE INSTRUCTIONS FOR USE, CONTACT THE COMPANY'S TECHNICAL SUPPORT.

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