Potenza FineX DISC

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PHS DO BRASIL LTDA

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Versão: 00

Potenza FineX DISC

INSTRUCTIONS FOR PROFESSIONAL USE

PHS DO BRASIL LTDA

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: 8.11992.7 (X401M5619X7L)

Notification at ANVISA number: 81199270034

Product Commercial Name: Potenza FineX Disc

Alternative Commercial Name: Abrasive Discs

Technical Product Name: Abrasive discs for finishing and polishing.

READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.

KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.

Technical Support

0800-000-3780

www.phsgroup.com.br

BASIC COMPOSITION:

The POTENZA FINEX DISC is composed of polyester film, aluminum oxide abrasive grains, an adhesive layer, and a metal eyelet made of a copper-zinc alloy, with a surface coating of tin and epoxy resin.

The POTENZA FINEX DISC KIT includes a mandrel made of AISI 304 stainless steel.

PRODUCT DESCRIPTION:

The POTENZA FINEX DISC is a system of dental abrasive discs intended for the reduction, finishing, and pre-polishing of restorative materials, such as composite resins, cements, and other restorative materials compatible with its intended use.

The discs are made of polyester film coated with aluminum oxide abrasive grains and feature a central metal eyelet for attachment to the POTENZA FINEX MANDREL, allowing their use with low-speed dental contra-angles.

The product is available in four grit sizes — Extra Coarse (80 µm), Coarse (60 µm), Medium (40 µm), and Fine (20 µm) — identified by colors, enabling the sequential performance of the reduction, finishing, and pre-polishing stages of the restorative surface.

The discs are for single use and must be discarded after the procedure. The POTENZA FINEX MANDREL is a reusable accessory, supplied in product kits or sold separately, and must be cleaned and sterilized before each use, according to its Instructions for Use.

PRODUCT' INDICATION:

The POTENZA FINEX DISC abrasive disc system is intended for the reduction, finishing, and polishing of composite resin dental restorations, removing surface irregularities, adjusting cervical margins, and providing progressive smoothing of the texture until pre-polishing is achieved. The product is designed for professional use by dentists, attached to a standard mandrel for a low-speed handpiece.

PRECAUTIONS:

- Before each use, check that the POTENZA FINEX DISC and the mandrel do not show signs of damage, deformation, or improper fixation. Do not use damaged components.
- Make sure that the disc is completely seated on the mandrel before starting rotation.
- Use only dental contra-angles in good working condition and compatible with the mandrel.
- Observe the maximum recommended rotational speed for each grit size. Use at higher speeds may cause premature wear, disc detachment, or reduced performance.
- Apply light and continuous pressure during the procedure. Excessive force may cause heating of the restoration, unnecessary wear of the disc, and reduction of its service life.
- Continuous irrigation during use is recommended to assist in dissipating the heat generated by friction and to facilitate the removal of debris produced during abrasion.

- Use the grit sizes in the recommended sequence, progressing from the coarsest to the finest, in order to obtain a better surface finish.
- If the disc shows vibration, instability, detachment, or any abnormal behavior during use, stop the procedure immediately and replace the product.
- Avoid unnecessary contact of the rotating disc with soft tissues, adjacent restorations, or dental structures that are not being treated.
- The product must be used only by dentists or legally qualified professionals trained in restorative dental procedures.
- Do not use the product if the packaging is damaged or has been opened.
- Keep the product in its original packaging until the time of use, preserving its characteristics and preventing contamination.

CONTRAINDICATIONS:

- The product must not be used in patients with known hypersensitivity or allergy to any of its components.
- The product must not be used directly on soft tissues, gingiva, oral mucosa, or skin.
- The product is contraindicated for any use other than the intended purpose established by the manufacturer, as described in these Instructions for Use.

POSSIBLE SIDE EFFECTS:

- Transient tooth sensitivity after the procedure.
- Temporary irritation of the gingiva or oral mucosa in the event of accidental contact with the rotating disc.
- Heating of the restoration or dental structure in the event of excessive rotational speed, excessive pressure, or insufficient irrigation, which may cause transient sensitivity.
- Mild discomfort during or after the procedure.

WARNING:

- Product for professional use only.
- The abrasive discs are for single use. Do not reuse. Reuse may result in loss of performance and cross-contamination.
- The POTENZA FINEX MANDREL is a reusable accessory and must be cleaned and sterilized before each use, according to its Instructions for Use. Failure to properly clean and sterilize the mandrel may result in cross-contamination.
- Do not exceed the maximum recommended rotational speed for each grit size. Use at higher speeds may result in excessive heating of the restoration or dental structure, as well as increase the risk of device failure.

- The use of proper isolation of the operative field is recommended to reduce the risk of accidental swallowing or aspiration of the disc or its fragments.
- Wear eye protection for both the professional and the patient during the procedure to protect against projected particles.
- Wear a protective mask to reduce the risk of inhalation of particles generated during abrasion.
- The use of a dental barrier is recommended, when indicated by the clinical technique, to help prevent accidental ingestion or aspiration of materials.
- Rotary instruments must be handled carefully and with control, avoiding abrupt movements.
- Inadvertent contact of the rotating disc with soft tissues (gingiva, oral mucosa, or skin) may result in local injury.
- Immediately discontinue use if the disc shows deformation, fracture, abrasive detachment, loosening of the eyelet, or any other sign of damage.

INSTRUCTIONS FOR USE:

The use of the POTENZA FINEX MANDREL, supplied in the kit or sold separately, is recommended, as it provides proper fixation and performance of the POTENZA FINEX DISC.

Step 1: Position the POTENZA FINEX DISC on the mandrel with the abrasive surface facing the surface to be treated. The abrasive surface corresponds to the side where the eyelet has a smaller thickness, while the side where the eyelet has a greater thickness corresponds to the back of the disc.

Step 2: Ensure that the disc is correctly seated and firmly attached to the mandrel before use.

Step 3: Insert the shank of the POTENZA FINEX MANDREL into a compatible dental contra-angle. Before starting the procedure, perform a brief functional test to verify the stability of the assembly and the proper rotation of the disc.

NOTE: The POTENZA FINEX DISC is intended for the reduction, finishing, and pre-polishing of dental restorations made of composite resins, cements, and other restorative materials, and should be used sequentially, starting with the coarsest grit and progressing to the finer grits according to the clinical need. Final polishing of the surface should be performed using specific polishing systems, when indicated. Use continuous water irrigation throughout the procedure to minimize heating of the treated surface, preserve dental tissues, and facilitate the removal of debris generated during abrasion.

Step 4: Apply the disc using light pressure and continuous, intermittent movements, avoiding excessive pressure on the treated surface.

Step 5: The following grit sequence is recommended:

- ✓ Extra Coarse 80 µm – (Purple): Initial reduction and removal of more pronounced excess material.
- ✓ Coarse 60 µm – (Yellow): Initial surface leveling.
- ✓ Medium 40 µm – (Blue): Surface finishing and refinement.
- ✓ Fine 20 µm – (Pink): Final finishing and surface preparation for the polishing stage.

IMPORTANT: Extra Coarse and Coarse grit discs must be used at a maximum speed of 12,000 rpm. Medium and Fine grit discs must be used at a maximum speed of 18,000 rpm.

Step 6: The procedure should be interrupted between the use of different grit sizes to rinse and dry the treated surface, allowing evaluation of the finish obtained before proceeding to the next stage.

Step 7: After completing the procedure, remove the POTENZA FINEX DISC from the mandrel and dispose of the abrasive disc in accordance with the applicable legislation.

CONSERVATION AND STORAGE:

- Store the POTENZA FINEX DISC in a clean, dry, well-ventilated place, protected from moisture, dust, direct sunlight, and corrosive chemical agents.
- Store the product at a temperature between 32°F and 104°F (0°C and 40°C).
- Keep the product in its original packaging until the time of use, preserving its physical integrity and preventing contamination.
- Avoid bending, creasing, compressing, or deforming the discs, as well as exposing them to impacts or contact with objects that may damage the abrasive surface or the metal eyelet.

TRANSPORT CONDITIONS:

- The product must be transported in its original packaging to preserve its physical integrity and prevent deformation, mechanical impacts, punctures, contamination, and exposure to corrosive chemical agents.
- During transportation, the POTENZA FINEX DISC must be maintained at a temperature between 32°F and 140°F (0°C and 60°C), protected from excessive moisture and direct exposure to adverse weather conditions.
- When applicable, stacking shall comply with the packaging strength limits, avoiding excessive compression that could deform the discs or damage the packaging.

DATE OF MANUFACTURE: SEE PACKAGE.

SHELF LIFE: SEE PACKAGE.

FOR DENTAL OFFICE USE ONLY. FOR PROFESSIONAL USE ONLY.

ENSURE THAT THESE INSTRUCTIONS FOR USE CORRESPOND TO THE PRODUCT BEING USED.

TO OBTAIN A PRINTED COPY OF THESE INSTRUCTIONS FOR USE, PLEASE CONTACT THE COMPANY'S TECHNICAL SUPPORT.

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