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## MANDRIL

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**PHS DO BRASIL LTDA**

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**Versão: 01**

# Potenza FineX MANDRIL

INSTRUCTIONS FOR PROFESSIONAL USE

**PHS DO BRASIL LTDA**

Guilherme Boldt Street, 210, District: Zona Industrial Norte, City: Joinville,

State: Santa Catarina, Country: Brazil - Zip Code: 89239-260

CNPJ: 17.910.389/0001-05 - BRAZILIAN INDUSTRY

E-MAIL: contato@phsdo brasil.com.br

Health Ministry's Operating Licence: 8.11992.7 (X401M5619X7L)

Notification at ANVISA number: 81199270033

**Product Commercial Name:** Potenza FineX Mandril

**Alternative Commercial Name:** Mandril

**Technical Product Name:** Mandril

**READ ALL THE INFORMATION IN THIS INSTRUCTION FOR USE CAREFULLY BEFORE APPLYING THE PRODUCT.**

**KEEP THIS INSTRUCTION UNTIL THERE IS NO FURTHER INTERACTION OF THE PRODUCT WITH YOUR PATIENT.**

**Technical Support**

**0800-000-3780**

[www.phsgroup.com.br](http://www.phsgroup.com.br)

**BASIC COMPOSITION:**

AISI 303 Stainless Steel.

**PRODUCT DESCRIPTION:**

The POTENZA FINEX MANDRIL is a reusable dental medical device, manufactured from AISI 303 stainless steel, designed for use with low-speed dental contra-angles. The product features a fixation system intended for coupling the POTENZA FINEX DISC, providing its retention during use.

The POTENZA FINEX MANDRIL acts exclusively as a mechanical element for supporting, fixing, and transmitting rotational movement to the POTENZA FINEX DISC, having no therapeutic purpose of its own and no direct clinical action on the patient.

The product is reusable and must undergo cleaning and sterilization processes before each reuse, as described in these Instructions for Use.

**PRODUCT' INDICATION:**

The POTENZA FINEX MANDRIL is indicated for the fixation and operation of the POTENZA FINEX DISC during dental procedures for reduction, finishing, and pre-polishing of restorative materials when used with low-speed dental contra-angles. The product is intended for professional use by dentists.

**PRECAUTIONS:**

- Before each use, verify that the POTENZA FINEX MANDRIL does not show signs of wear, deformation, corrosion, cracks, or any other damage that may compromise its use.
- Ensure that the POTENZA FINEX DISC is correctly seated and firmly attached to the mandrel before starting the procedure.
- Use only compatible dental contra-angles that are in good working condition.
- Before each use, perform a brief functional test to verify the stability of the assembly and the proper rotation of the disc.
- Do not use the mandrel if it exhibits vibration, instability, or difficulty in securing the disc.
- After each use, the mandrel must undergo cleaning, drying, inspection, and sterilization before being reused

**CONTRAINDICATIONS:**

- The product is contraindicated for any use other than the intended purpose established by the manufacturer, as described in these Instructions for Use.
- Do not use the product if it shows deformation, corrosion, cracks, or any other damage that may compromise its safety or performance.

**POSSIBLE SIDE EFFECTS:**

- No side effects are known when the product is used in accordance with these Instructions for Use.

**WARNING:**

- Product for professional use only.
- The POTENZA FINEX MANDRIL is reusable and must be cleaned and sterilized before each use, as described in these Instructions for use.
- Do not use the product if it shows signs of wear, deformation, corrosion, cracks, or any other damage that may compromise its integrity.
- Do not use the product unless the POTENZA FINEX DISC is correctly attached to the mandrel.
- Rotary instruments must be handled carefully and with control, avoiding abrupt movements.
- Inadvertent contact of the rotating mandrel with soft tissues (gingiva, oral mucosa, or skin) may result in local injury.

**INSTRUCTIONS FOR USE:**

**Step 1:** Before use, verify that the POTENZA FINEX MANDRIL shows no signs of wear, deformation, corrosion, cracks, or any other damage that may compromise its use.

**Step 2:** Attach the abrasive disc to the mandrel fixation system, ensuring that it is correctly seated and firmly secured.

**Step 3:** Insert the shank of the POTENZA FINEX MANDRIL into a compatible dental contra-angle and ensure that it is properly secured.

**Step 4:** Before starting the procedure, perform a brief functional test to verify the stability of the assembly and the proper rotation of the abrasive disc.

**Step 5:** Perform the clinical procedure according to the Instructions for Use of the abrasive disc being used.

**Step 6:** After completing the procedure, remove the abrasive disc from the mandrel and dispose of it according to the disc manufacturer's Instructions for Use.

**Step 7:** Subject the POTENZA FINEX MANDRIL to the cleaning, drying, inspection, sterilization packaging, and sterilization processes before reuse, as described in these Instructions for Use.

**Note:** The use of the POTENZA FINEX DISC abrasive disc is recommended to ensure compatibility, proper fixation, and the performance intended by the manufacturer.

## **CLEANING, DISINFECTION, AND STERILIZATION**

The Potenza FineX Mandril is supplied in a non-sterile condition and must be cleaned and sterilized between successive uses.

### **Cleaning**

Immediately after use, remove visible residues under running water.

Wash the instrument using an enzymatic detergent according to the detergent manufacturer's instructions, using a soft-bristle brush when necessary to remove adhered residues.

Rinse thoroughly under running water until all detergent residues have been completely removed.

Dry the instrument completely using a clean, lint-free cloth or oil-free, moisture-free compressed air.

### **Sterilization**

After drying, place the mandrel in suitable steam sterilization packaging, properly sealed.

Sterilize in a steam autoclave at 250°F (121°C) for 20 to 30 minutes or, alternatively, at 273°F (134°C) for 5 to 10 minutes, according to the protocol validated by the sterilization service.

At the end of the cycle, allow the packaging to dry before removing it from the autoclave.

Store the sterilized product in a clean, dry place, protected from moisture and direct sunlight, preserving the integrity of the packaging until the time of use.

## **CONSERVATION AND STORAGE:**

- Store the POTENZA FINEX MANDRIL in a clean, dry, well-ventilated place, protected from moisture, dust, and corrosive chemical agents.
- Store the product at a temperature between 32°F and 104°F (0°C and 40°C).
- After cleaning, drying, and sterilization, keep the mandrel packaged in suitable steam sterilization packaging, preserving the integrity of the sterilization packaging and the sterility of the product until the time of use.

## **TRANSPORT CONDITIONS:**

- The product must be transported in its original packaging to preserve its physical integrity and prevent deformation, mechanical impacts, punctures, contamination, and exposure to corrosive chemical agents.
- During transportation, the POTENZA FINEX MANDRIL must be maintained at a temperature between 32°F and 140°F (0°C and 60°C), protected from excessive moisture and direct exposure to adverse weather conditions.

- When applicable, stacking shall comply with the packaging strength limits, avoiding excessive compression that could damage the product or its packaging.

**DATE OF MANUFACTURE: SEE PACKAGE.**

**SHELF LIFE: SEE PACKAGE.**

**FOR DENTAL OFFICE USE ONLY. FOR USE BY DENTAL PROFESSIONALS ONLY.**

**ENSURE THAT THESE INSTRUCTIONS FOR USE CORRESPOND TO THE PRODUCT BEING USED.**

**TO OBTAIN A PRINTED COPY OF THESE INSTRUCTIONS FOR USE, PLEASE CONTACT THE COMPANY'S TECHNICAL SUPPORT.**

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